



Environment, soil, geology, marine science



Organic chemistry, pharmacy



Food analysis & security



Petrol chemistry, coal, energy




Quality control


Carbon, Nitrogen and Hydrogen in coal


Instrument: ECS 8020


Mode: CHNS

Pretreatments: none

 Via Milano, 15/A - 20060 Bussero (MI), Italy

 Phone: +39 02 950 34 69

 www.nctechnologies.it

 info@nctechnologies.it



NC TECHNOLOGIES

Innovative Elemental μ -Analysis

Carbon, Nitrogen and Hydrogen in coal



The chemical and physical characterization of coals is important for evaluating their use as fuels for the production of electricity.

Based on their degree of aging, the coals have different percentages of Carbon and are classified in:

- peat (C from 55 to 65%)
- lignite (C from 65 to 75%)
- bituminous carbon (C from 75 to 90%)
- anthracites (C > 90%)


The characterization of coal occurs mainly through different analytical methods like thermal analysis, ashes analysis and elemental analysis.

ECS8020 can perform reliable analyses in a very short time; research and applications let us to optimize the instrument and the consumables for the perfect coal combustion and conversion to detectable gases.

Parameter	Carbon	Nitrogen	Hydrogen
Average	55.96	1.35	4.62
Standard deviation	0.90	0.02	0.07
Average accuracy	1.60	1.64	1.45

All reported values unit: %

- ✓ Configuration: CHNS
- ✓ Furnaces: no. 2
- ✓ Sampler: Pneumatic
- ✓ Chemical standard: Acetanilide



To send your samples for free demonstration analyses:
info@nctechnologies.it

For analytical and technical questions:
customerservice@nctechnologies.it

