

Environment, soil, geology, marine science

Organic chemistry, pharmacy

Food analysis & security

Petrol chemistry, coal, energy

**Quality control** 

## **Carbon and Nitrogen**

## in soils

Instrument: ECS 8024 Soil Special Mode: CN Pretreatments: none

Via Milano,15/A - 20060 Bussero (MI), Italy Phone: +39 02 950 34 69 www.nctechnologies.it info@nctechnologies.it



## Carbon and Nitrogen in soils



Soil Sample	Carbon	Nitrogen
Cultivated soil 1	$0.632 \pm 0.012$	0.065±0.003
Cultivated soil 2	$1.686 \pm 0.044$	0.122±0.075
Uncoltivated soil	3.411±0.013	0.474±0.008
Average accuracy	1.629	2.348
All reported values unit: %		


Carbon and Nitrogen determination in soils are fundamental data for soil chemistry characterization, environmental monitoring, agronomy and cultures.

Carbon and Nitrogen are diffused elements in soils but are normally reported into low concentrations. Considering the cultivated soils or poor natural soils, concentrations of these elements are usually lower than 2% and 0.2% for carbon and nitrogen respectively (typical C:N ratio reported in cultivated soils is 10).

Low elements concentration is also accompanied by the presence of a pool of different minerals, metals, organic compounds that may cause interferences or problems in detection.

ECS 8024 Soil Special is the optimized analyzer for CHNS-O determination in any type of soil. Consumables are set and tested for the perfect oxidation and reduction of sought elements.

Thanks to the optimized conditions and to the equilibrium in consumables and chemical choice, very high analysis accuracy and precision are ensured. The repeatability allow the user to save time, chemicals, products and thus money.



To send your samples for free demonstration analyses: info@nctechnologies.it For analytical and technical questions: customerservice@nctechnologies.it



