



Agri - Food



Chemicals



Geo - Environment



Energy

Application Note

Biogenic CO₂ in atmospheric air

Instrument: 8070 AIR CO₂ – 14C SCAR

Mode: C



Via Milano, 15/A - 20060 Bussero (MI), Italy



Phone: +39 02 950 34 69



www.nctechnologies.it



info@nctechnologies.it



NC TECHNOLOGIES

Innovative Elemental Analysis

Application Note – Biogenic CO₂ in atmospheric air

The Greenhouse Gas (GHG)

Scientific evidence and economic statistics indicate that increases in the abundance of atmospheric greenhouse gases (GHGs) since the Industrial Revolution are the result of human activities, either directly or indirectly. Furthermore, the convergence of extensive scientific evidence indicates that increases in long-lived, well-mixed GHG, particularly CO₂, CH₄, nitrous oxide (N₂O), are the main cause for increases in global temperature over the industrial period.


Increasing of CO₂

There are both natural and human sources of carbon dioxide emissions. Natural sources include decomposition, ocean release and respiration. Human sources come from activities like cement and steel production, deforestation, waste incineration plants as well as the burning of fossil fuels like coal, oil and natural gas.

Due to human activities, the atmospheric concentration of carbon dioxide has been rising extensively since the Industrial Revolution and has now reached dangerous levels not seen in the last 3 million years.

Human sources of carbon dioxide emissions are much smaller than natural emissions but they have upset the natural balance that existed for many thousands of years before the influence of humans.

This is because natural sinks remove around the same quantity of carbon dioxide from the atmosphere than are produced by natural sources. This had kept carbon dioxide levels balanced and in a safe range. But human sources of emissions have upset the natural balance by adding extra carbon dioxide to the atmosphere without removing any.



Annual mean carbon dioxide growth rates based on globally averaged marine surface data.

Data by The Global Monitoring Laboratory (GML) of the National Oceanic and Atmospheric Administration (NOAA)

Application Note – Biogenic CO₂ in atmospheric air


Fossil CO₂ emission sources

Animals – including humans – breathe out CO₂ in the process of respiration. When plants and animals die, their stored carbon is also released as CO₂. Natural processes such as respiration and decay, forest fires and volcanic eruptions add an additional 190.2 billion tonnes of CO₂ to the atmosphere per year.

Excess carbon dioxide emitted by human activity is impacting the natural carbon cycle in the environment that has occurred for millions of years.

Humans generate CO₂ when burning fossil fuels such as gas, petrol, oil, and coal. This adds an additional 9.1 billion tonnes of CO₂ to the atmosphere each year. Plants and soils take up 2.8 billion tonnes of this extra carbon, while the oceans take up 2.2 billion tonnes.

The remaining 4.1 billion tonnes of CO₂ stays in the air, increasing the atmospheric concentration of CO₂.



CO2 in air worldwide. Data by NASA

About 90 per cent of the world's carbon emissions comes from the burning of fossil fuels – mainly for electricity, heat and transport.

In 2022, most of the world's fossil fuel carbon emissions came from coal (40 per cent), oil (32 per cent), natural gas (21 per cent), cement (5 per cent) and flaring and other smaller sources (2 per cent).

Just four regions accounted for about two-thirds of global fossil-fuel carbon emissions in 2021: China (31 per cent), the USA (14 per cent), the EU27 (7 per cent), and India (7 per cent).

Application Note – Biogenic CO₂ in atmospheric air

Biogenic CO₂

When we talk about carbon dioxide (CO₂), it's easy to think of it as a single problem: a greenhouse gas driving climate change. But not all CO₂ is created equal. Its source, biogenic or fossil, makes a huge difference in how it impacts the environment, our climate targets, and the sustainability of industrial processes.

Biogenic CO₂ comes from natural, short-term carbon cycles. It's produced when organic materials, like plants, crops, or biomass, decompose or are combusted.

Biogenic CO₂ in atmospheric air refers to carbon dioxide released by biological sources, such as plant respiration, decomposition, and organic combustion. Because this CO₂ was recently sequestered from the air via photosynthesis, it is considered part of the natural, short-term carbon cycle and does not add new net carbon to the atmosphere.


The distinction between biological and fossil CO₂ is critical for climate accounting and mitigation strategies:

The Natural Carbon Cycle: plants and trees absorb atmospheric CO₂ through photosynthesis to grow. When these plants decompose, are consumed, or are burned as biofuels, they release that exact same CO₂ back into the atmosphere. This creates a continuous, circular biological loop rather than introducing new carbon into the environment.

Carbon Neutrality: because it relies on recently captured carbon, biogenic CO₂ is often treated as net-zero under the Greenhouse Gas Protocol and is exempt from many industrial trading frameworks like the European Union Emissions Trading System.

Distinguishing Sources: scientists use radiocarbon (¹⁴C) dating to measure atmospheric samples. Fossil fuels have been trapped underground for millions of years and contain no radiocarbon whereas biogenic CO₂ reflects current atmospheric radiocarbon levels, allowing researchers to accurately separate the two types of emissions

Applications: biogenic CO₂ is increasingly captured from organic waste and biomass processing, then combined with green hydrogen to synthesize renewable e-fuels and materials, offering a zero-greenhouse-gas-intensity alternative to fossil sources



Application Note – Biogenic CO₂ in atmospheric air

How to measure biogenic fraction of CO₂

When biomass is used as a fuel, carbon stored in renewable organic material is released into the atmosphere. It is captured back into the natural carbon cycle as the biomass regrows.

Burning a mixture of biomass and fossil fuel can significantly reduce fossil fuel consumption and lower the total carbon footprint.

When the two fuels are burned together, an accurate breakdown of what fraction of the emissions are biogenic CO₂ can be achieved using the process of radiocarbon dating.

This is possible because biogenic CO₂ emissions contain carbon-14 isotopes and fossil fuel CO₂ emissions do not.

Practical implications of studying the biogenic fraction of CO₂

The determination of the biogenic fraction of atmospheric CO₂ is of fundamental importance for assessing air quality and monitoring its temporal evolution. However, although such investigations have traditionally been confined to the research community, they are not the only applications of scientific and practical relevance.

Indeed, the determination of the biogenic CO₂ fraction also has significant implications for a range of real-world applications beyond air quality studies.

Companies with large amounts of emissions arising from partly biogenic materials can greatly benefit by differentiating the fraction that is attributable to biogenic versus fossil CO₂.



Application Note – Biogenic CO₂ in atmospheric air

The benefits include:

- capitalising on potential surplus allowances gained from more precise quantification of fossil CO₂;
- subtracting biogenic emissions results in lower annual emissions for reporting and potentially generates surplus allowances for trading;
- separate determination of biogenic and fossil CO₂ emissions helps organisations quantify their true GHG inventory; and
- accurate accounting of fossil CO₂ emissions means that organisations know the correct amount of credits needed to offset their carbon footprint over a certain period.

The biogenic radiocarbon testing

Due to the many uncertainty factors associated with mixed-fuel sampling and analysis, radiocarbon dating is recognised as a much more accurate and cost-effective way to determine biogenic CO₂.


It also:

- ✓ allows for fast and easy sample extraction or installation of monitoring equipment
- ✓ turns around results quickly
- ✓ accepted under EU Emissions Trading Scheme (ETS), California Air Resources Board (CARB), the US/ Canadian Western Climate Initiative (WCI) and the Australian National Greenhouse and Energy Reporting (NGER) Program.

Anthropogenic CO₂ emissions from fossil fuels play an important role in the environmental pollution.

Measurements of the radioactive carbon isotope ¹⁴C in CO₂ are the only quantitative observational means to decipher the fossil fuel CO₂ component in the atmosphere.

This is because CO₂ from burning of fossil fuels, due to the large age of coal, oil and gas of several hundred million years, is free of ¹⁴C. Adding fossil fuel CO₂ to the atmosphere, therefore, not only leads to an increase of its mixing ratio but also to a decrease of the ¹⁴C/C ratio in atmospheric CO₂.




Application Note – Biogenic CO₂ in atmospheric air

A new and completely automatic ¹⁴C analysis method

The **¹⁴C SCAR** (Saturated-absorption **CA**vity **R**ing-down) spectrometer is a new and innovative spectrometer for the ¹⁴C measure.

Developed by a team of the Italian National Research Council, It is able to measure the mole fraction of radiocarbon in any sample using the SCAR technique, which improves the CRD (cavity ring down) limits by more than an order of magnitude.



Schematic of the optical and optoelectronic setup of the radiocarbon detection part of C14-SCAR.

In the opto-electronic setup, two Quantum Cascade Lasers (QCLs) are used as laser sources.

The QCL2 is frequency-locked to the transition of the N₂O molecule, the laser thus locked to the molecular line, via direct spectroscopy, is used as an absolute reference for the frequency scanning.

The QCL1 is tuned, thanks to the Pound-Drever-Hall locking (PDH) technique, to the resonant cavity (Fabry-Perot cavity) and the laser used to perform the SCAR spectroscopy of the transition of ¹⁴C¹⁶O₂.

A further frequency-locking occurs between the beat note of the two QCLs, fixing the frequency of the QCL1 to the reference one (QCL2).

Patent [WO2014170828A1](#), "Apparatus and method for measuring the concentration of trace gases by SCAR spectroscopy".

Patent [WO2017055606A1](#), "Method for measuring the concentration of trace gases by SCAR spectroscopy"

Patent [WO2016067241A1](#), "Low-noise current source"



Compliance to the UNI EN 16640

Bio-based products –Bio-based carbon content –

Determination of the bio-based carbon content using the radiocarbon method



Application Note – Biogenic CO₂ in atmospheric air

14C isotope analysis method comparison

The detection of 14C is really challenging , due to its extremely low concentration in nature: the required sensitivity is of 1 part in 10¹⁵

Nowadays, there are three main analytical methods capable of achieving the required sensitivity:

- ❖ Accelerator Mass Spectrometry (AMS)
- ❖ Liquid Scintillation Counting (LSC).
- ❖ Saturated-absorption CAvity Ring-down (14C SCAR)

Below is a brief summary comparison table

	14C - SCAR	AMS	LSC
Size	 	  	 
Weight	 	     	 
Minimum sample needs	 	 	     
Precision measurement	 HIGH	 HIGH	 LOW
Dynamic range	 VERY HIGH	 LOW	 HIGH
Sample handling and preparation	 Fully automatic	 automatic	 manual

Application Note – Biogenic CO₂ in atmospheric air


Not only 14C SCAR

For this specific application, the 14C SCAR detector was coupled to a sampling unit capable of collecting ambient air, removing all gases except CO₂, and concentrating the CO₂ prior to analysis.

8070 Air CO₂ is a concentrator and purifier of CO₂ contained in atmospheric air. Therefore, through its use, it is possible to obtain purified CO₂ and introduce it into mass detectors to determine the isotopic ratio $\delta^{13}\text{C}$ or quantify the content of 14C, a radioactive isotope of carbon.



8070 Air CO₂



14C - SCAR


This unit is based on the production of a high volume of pure CO₂ according to the principle of adsorption/desorption, with an innovative purification line that allows the elimination of water, VOC and NO_x, leaving only pure CO₂.

With the new C-Quantum CO₂ absorption system, a large amount of carbon dioxide can be treated with an automatic regeneration system, achieving more precise results and higher performance than other systems. The 8070 Air CO₂ is particularly suitable for connection to other units to determine the isotope ratios of stable carbon isotopes or radiogenic 14C in a fully automatic way

Application Note – Biogenic CO₂ in atmospheric air


The integration of the 8070 Air CO₂ system with the 14C SCAR detector enabled fully automated atmospheric air sampling, eliminating the need for cryogenic processing.

The graph below shows atmospheric air samples collected at hourly intervals over a total sampling period of approximately 12 hours. It can be observed that the increase in atmospheric CO₂ concentration is accompanied by a corresponding decrease in the radiocarbon (14C) content, consistent with the contribution of fossil fuel emissions.



The relationship between carbon dioxide (CO₂) and radiocarbon (14C) is essential for tracking pollution. Radiocarbon (14C) is absent from fossil fuels (coal and petroleum). When these fuels are burned, they release CO₂ that is devoid of radiocarbon, reducing the proportion of 14C in the atmosphere and providing evidence that human activities are responsible for climate change.


The graph below shows atmospheric air samples collected at hourly intervals over a total sampling period of approximately 24 hours.




Application Note – Biogenic CO₂ in atmospheric air

The graph below presents the results of an analysis of the exhaust emissions from a typical automobile.

A low concentration of carbon-14 (¹⁴C) indicates that the fuel used is predominantly of fossil origin.



Atmospheric reference value: 98.5 pMC



Biobased Carbon Content: 7%

Conclusions

The 8070 Air CO₂ coupled with the 14C SCAR is a highly promising and innovative system for analyzing biogenic content in atmospheric CO₂.

Its:

- ❖ performance
- ❖ automatism
- ❖ compact size
- ❖ ease of use
- ❖ no real sample preparation

make it particularly attractive for a wide range of commercial, research, and development applications.